


A Magnetic Resonance Imaging (MRI) System is designed to “see” inside a given object or tissue sample (living or dead), and produce an image that represents, as closely as practical, how the object or sample really looks on the inside. This image can be used then to diagnose physical and chemical abnormalities in the sample – or indeed declare the sample to be free of abnormalities (i.e. normal). MRI techniques are now commonly used for medical diagnoses, materials research, and newer industrial applications (test food quality, test timber quality, rock permeability, etc.). In order to do this, an MRI machine uses a combination of one large stationary magnet, three directional and precisely controlled electro-magnetic coils (gradient coils), an RF pulse generator, RF detector, and the properties of elemental nuclear magnetism.

Inside the nucleus of every atom, individual protons and neutrons spin about an axis. This property, called spin angular momentum, is the basis of nuclear magnetism. Since atomic nuclei have charge, this spinning motion produces a magnetic moment along the spin axis. In most nuclei, the particles are paired so that the net magnetic properties cancel within each atom. However, if the number of protons or neutrons is odd, complete cancellation is not possible. Nuclei with an unpaired proton or neutron such as hydrogen-1, carbon-13, and sodium-23, among others, exhibit a net magnetic effect. The relative strength of this magnetic moment is a unique property for each element, and therefore determines its suitability for magnetic resonance (MR) energy absorption and detection. **The hydrogen (^1H) nucleus, which is highly abundant in biological systems, has the strongest magnetic moment.** In MRI for human tissue then, the system takes advantage of the specific nuclear magnetic properties of our hydrogen atoms, and how those nuclear magnetic properties interact with large, external, controlled magnetic fields.

If you look at this basic MRI system block diagram, you can get an idea of how the various blocks work together to produce a usable image. The huge static magnet, usually providing a magnetic field measuring either 1.5T (Tesla) or 3.0T (Tesla), is shown in blue. It is usually enclosed in liquid helium and the enclosed unit often weighs between 9 and 12 tons. This large superconducting electromagnet provides the strong, stable magnetic field that lines-up the nuclear magnetic spins of the H-atoms inside the body.



An Introduction to MRI (cont.)

The real MRI magic is in the RF sections (lime green) and the magnetic field Gradient sections (red) which work together respectively to: 1. Set-up the magnetic resonance conditions within the “slice” of the sample under test at a given moment; and, 2. Encode those specific resonating atoms with three-dimensional “markers” to pinpoint the location of the tissue being analyzed (in real time). The transmit/receive sequence (roughly stated here) starts with a particular individual RF coil transmitting a measured pulse of RF energy at the specific 1H Larmor frequency (F_0). For a 1.5T machine, F_0 is 63.864 MHz; For a 3.0T machine, F_0 equals 127.728 MHz. At the end of the Transmit Pulse, the 1H atoms are resonating in their new axis. Next the Gradient coils are pulsed to encode the X-Y-Z coordinates into the resonance, so that shortly thereafter, a receive RF Coil can detect the tiny RF signal (frequency and encoding) being released as that particular group of atoms head back toward their previous state. In this way, detailed information is gathered (in an iterative fashion) about the precise chemical make-up of each small area of the tissue being resonated, sample section by sample section.

The entire system is controlled by a sophisticated computer running specialized software (including different programs for each type of scan - full-body, knee, ankle, brain, heart, etc.). The software controls the precise sequencing for each dynamically controlled sub-system (all RF sections, and the X, Y, and Z gradient coils), and collects/analyzes the large amount of received data to product an image. The computer even tunes the RF coils, in real time at the very beginning of each MRI session, by impedance-matching the coils to the sample under test (tuning the RF Coils for maximum Signal-to-Noise Ratio).

Contact your local sales representative or learn more about Richardson RFPD online at www.richardsonrfpd.com.

Your Global Source for
RF, Wireless & Energy Technologies

